

microlife®

microlife®

WatchBP Vital System



Microlife US Office
1617 Gulf to Bay Blvd., 2nd Floor
Clearwater, FL 33755
www.microlife-us.com

1-866-950-1404

Contact mail: contact@microlifeusa.com

Benefits

- Reduce the Time Required to Collect Vital Signs
- Auto-Populate Vital Data with Date and Time Stamp
- Increase the Efficiency & Effectiveness of collecting daily measurements
- Reduce Manual Key-in Error
- HIPPA Compliant
- Fully Integrated with PointClickCare®



PointClickCare®





Ease of Use

Learn to use the WatchBP Vital System within thirty minutes or less. Collected vitals will be sent directly to the Microlife WatchBP Vital App, once the collection of measurements has been completed.

Time Saving for Caregivers

Once the process has been completed, the collected vitals are automatically transferred to your PointClickCare® EHR platform eliminating the need for manual data entry. Decrease your facility's costs by increasing your staff's overall accuracy, efficiency and effectiveness. Microlife WatchBP Vital System allows a hassle-free setup.

Easy to Maintain


Microlife WatchBP Vital System connects with the Blood Pressure Monitor, Blood Glucose Meter, Forehead Thermometer, and Oximeter via Bluetooth Connectivity. Greatly reduce downtime for maintenance and repair.

Customize Setting

Caregivers are able to customize features within the Microlife WatchBP Vital System App to standardize their workflows



microlife®



OXY 560 BT

Fingertip Pulse Oximeter

Measurement Range:

- SpO2: 0% - 100%
- Pulse Rate: 30 - 250 bpm
- PI: 0% - 20%



FR1MZ1-B

Forehead Thermometer

- Response Time: 1-2 Seconds
- Measurement Range: 32.2 - 211.8 F
- Non-contact Measurement (Distance within 2 inches)

BP M1

Blood Pressure Monitor with AFIB Screening

Powered by WatchBP Algorithm

- Clinically Validated with Diabetic Patients
- End Stage Renal Disease (ESRD) Patients
- AFIB Screening

